



Send your data to your doctor!

It is essential for optimal diabetes treatment that the doctor knows your blood sugar values and also knows when and how much insulin you are injecting and which medicines you are taking.

>> It is often very difficult to read such information from written paper diaries. Values are often missing there or are no longer legible.

Many doctors therefore read their patients' measurement devices into their computers. However, this is often laborious and time-consuming during the consultation.

It would therefore be more efficient to send the values before you visit your doctor. This will help your doctor to save time and also to prepare his consultation with you better.

This is all very easy, secure and free with gluconet®, the innovative telemedicine solution for diabetes patients.

Questions about gluconet®

What is gluconet®?

gluconet® is an innovative telemedicine solution specifically for diabetics. You can send your values from your measurement device directly to your doctor from your home computer.

Why all this?

Many doctors read their patients' measurement devices into their computers. However, this is often laborious and time-consuming. It would therefore be more efficient to send the values before you visit your doctor. This will help your doctor to save time and also to prepare your consultation better.

Why is gluconet® so secure?

Apart from the encrypted measurement values, you only need to send a number (the gluconet® ID) for identification purposes. It is not necessary to enter any personal data, not even your e-mail address. Even if someone succeeded in hacking into the data transfer and cracking the encryption, it would still not be possible to discover your identity from this.

Which blood sugar measurement device can I use?

gluconet® supports all common blood sugar measurement devices; you can also send data from diary apps, insulin pumps or blood pressure measurement equipment. You can find an up-to-date list under www.gluconet.org

Is it really free?

gluconet® is entirely free for you as a patient. You are not entering into any obligation and you are not even required to give your e-mail address.

What do I need?

You only need a computer with an internet connection and also a suitable data cable for your measurement device.

What is behind this?

gluconet® is operated by mediaspects GmbH, the leading Europe-wide producer of diabetes management software. DIABASS® PRO software is used in 1000s of medical practices and clinics, and more than 100,000 patients use DIABASS® for their daily diaries.




Telemedicine for diabetes patients

Send your blood sugar values from your home to your doctor

- Simple
- Secure
- Universal
- Free of charge



gluconet® is that simple:



You will receive a gluconet® Pass from your doctor which contains a short introduction and your personal gluconet® ID – this is a code which your doctor uses to classify your data.



Enter your gluconet® ID in www.gluconet.org and download the free data transfer software. No other data will be requested, nor is any registration required.



With a few clicks you can read the values from all common blood sugar measurement devices, blood pressure measurement devices and also insulin pumps into your computer, which are then automatically forwarded to your doctor. It is not necessary to enter an e-mail address or any other personal data.

- >> It goes without saying that your data will nevertheless be secured with a high level of encryption (AES256).
- >> All this is entirely free of charge – all you need is the data transfer cable for your measurement device.
- >> After the data transfer has finished you will receive a PDF document where your values will be displayed in a clear format as a diary. You can then print out and save this document.



By the way: gluconet® can also be used with apps. Just send an e-mail with the data and your gluconet® ID in the subject line to doc@gluconet.org.

If your diary app is certified for gluconet® ID then the data will be sent just by pressing a button.



Receiving data is just as easy for your doctor: existing data can be downloaded automatically with one click, and then assigned to the patient – locally therefore on the doctor's PC.



gluconet® can be used with 200 blood sugar measurement devices!
You can find an up-to-date list on the internet under www.gluconet.org

